


# "ZENITH" 548/E Antibacterial

## PLIER STAPLER

### Materials:

- made entirely in painted metal with nickel-plated finishing
- lead free and cadmium-free enamel painting highly resisting wear and corrosion and with highly functional additive for releasing silver ions



### Technical / functional features:



- **Stapling:**



- **Staple type:**

6/4

6/6

- **"ZENITH" staple ref.:**

130/E - 130/Z  
130/LM - 130/LL

130/Z 6

- **Loading capacity (N° of staples):**

100

100


- **Stapling capacity (N° of sheets - 80 g/sqm):**

15

30

- **Max. stapling depth (for paper insertion):**

47 mm



- **Colour:**

DEWY GREEN

- **Packaging:**

box containing a stapler and  
a box of "ZENITH" staples ref. 130/E

- **Box dimensions:**

15,4 x 8,7 x 1,6 cm

- **Weight:**

240 g (net) / 285 g (gross, including a box of  
"ZENITH" staples ref. 130/E)

- **N° of items per carton:**

60 (packets of 10 pieces)


- **Carton dimensions:** 32,5 x 27 x 21 cm (AxBxC)

- **Carton weight:** 17,5 Kg

**Loading system:** rear loading, with extractable pusher

**Patented anti-jamming device.**

The plier directions are on the back of the packaging.



### Descriptive notes:

"ZENITH" stapler 548/E Antibacterial maintains all the technical and functional characteristics of the 548/E "ZENITH" stapler, leader in the office and used every day by millions of users who are satisfied with its quality, practicality, strength and endless capacity to function, even after years of continuous use. The special additive for releasing highly functional silver ions allows the formation of a protective film on the whole stapler, which is resistant to attacks by strains of **Gram-negative** and **Gram-positive bacteria** such as **Escherichia Coli** and **Staphilococcus Aureus**.

**Quality:** to maintain the high quality standard guaranteed by all "ZENITH" products, each plier stapler is tested piece by piece. **Life-guarantee for manufacturing defects.**

**In addition:** The very high antibacterial resistance, evaluated according to ISO 22196 - 2011, does not interfere with the performances of the stapler, which is therefore suitable to be used in environments that require specific prevention and resistance to bacterial proliferation.

"ZENITH" 548/E Antibacterial can also staple with the "ZENITH" staples ref. 130/Z6 (length of legs 6 mm) up to 30 sheets!